



## **nonStick** red

The new bipolar forceps from KLS Martin –  
with nonStick effect and a revolutionary  
ergonomic design



revolutionary  
ergonomic  
design

## nonStick red

The new bipolar forceps from KLS Martin –  
with nonStick effect and a revolutionary  
ergonomic design

„With the new KLS Martin bipolar nonStick red forceps, tissue sticking to the tips of the instrument is definitely a problem of the past.

Especially for operations involving extremely fine dissection (such as parotid and thyroid interventions or microvascular tissue transfer), the use of this new generation of bipolar forceps provides significant advantages due to their exceptional haptic quality: minimized risk of procedures, gentler operating technique and shorter operating times. The tips of the instrument need to be cleaned considerably less often during interventions and the cleaning process itself is much faster as well.

We simply wouldn't like to do without these forceps.“

Dr. Paul-Stefan Mauz (M.D.)  
Senior Consultant  
Tübingen University Hospital  
Department of Otorhinolaryngology



## nonStick red

The new bipolar forceps from KLS Martin – with nonStick effect and a revolutionary ergonomic design


### **The new generation of bipolar forceps**

Who hasn't encountered this everyday operating situation when working with conventional bipolar forceps? The tissue tends to stick to the tips of the instrument during coagulation and when withdrawing the instrument, the coagulated tissue is torn open again.

The result is renewed bleeding. Now, this unwelcome effect can be prevented by using KLS Martin nonStick red, the new and innovative generation of bipolar forceps – the coagulated tissue is no longer disrupted when opening the instrument.

What's more, the cumbersome and time-consuming process of cleaning the tips of the forceps during operations is largely eliminated at the same time. This means fewer interruptions of the surgical procedure and corresponding time savings.

## nonStick red



Coagulating produces so-called „hot spots“.




Thanks to the special silver alloy used for the tips of the forceps, the heat can be discharged quickly and effectively.




This prevents the tissue from sticking to the tips of the forceps.

## Standard forceps




Coagulating produces so-called „hot spots“.



As the heat cannot be discharged, overheating is the inevitable result.



Tissue tends to stick to the tips of the forceps. When opening the instrument, charred tissue is torn off in the process.




Thermography image of the nonStick red forceps

### Improved heat conductivity

A special feature does the trick – polished silver tips offering excellent heat conductivity. Thanks to this property, the tissue cannot adhere to the tips of the forceps.

The heat is discharged from the tips quickly and effectively. And what's more, the nonStick effect will last because the tips are not merely coated, but made of solid noble metal.



Thermography image of a standard bipolar forceps


### Heat conductivity standard forceps

The thermography image shown here demonstrates that in the case of a standard bipolar forceps, the heat becomes trapped at the tips of the forceps. The heat cannot be quickly and effectively transported away. The tips of the forceps become clogged with tissue as a result.

# nonStick red

## Variety of models

The new bipolar nonStick red forceps from KLS Martin are available in all usual sizes for the most diverse fields of use.



## The tips of the forceps

The tips are precision-made and available in different sizes for extremely accurate surgical use. By adjusting the pressure applied, the user can coagulate merely with the tips or by using the full surfaces.

The polished silver-metal tips are excellent heat conductors preventing tissue adhesion.




## marSIGHT

The special geometry of the forceps tips gives the surgeon an unobstructed view of the surgical site and the tissue to be coagulated. Moreover, the tips enable you to grasp tissue with utmost precision. Don't miss out on the experience!




The models featuring tips with marSIGHT geometry are clearly marked in this brochure.

## revolutionary ergonomic design




### Guiding mechanism

The guiding mechanism ensures the accurate, highly parallel closure of the pre-tensioned forceps and prevents „tip spread“ as well. This, in turn, prevents fatigue in the user’s hand and facilitates the dissecting process.



### Connection

The bipolar nonStick red forceps can be used with any commonly available HF generator. Special connecting cables are available for this purpose.



### marGRiP

The ergonomic handle surface provides for secure handling and prevents fatigue in the user’s hand.



### Shank geometry

Thanks to the view-optimized geometry of the shanks of the forceps, the surgeon’s hand and the cable connector are kept out of the surgical field. This is particularly advantageous when working under the operating microscope.

revolutionary  
ergonomic  
design

Bipolar forceps, KLS Martin nonStick red, straight



1/1

80-982-12-04  
12 cm/4 3/4"  
pointed, 0.3 mm


Example illustration

Tip width	Length	Item No.
<b>pointed</b>		
0.3 mm	12 cm/4 3/4"	80-982-12-04
0.3 mm	17 cm/6 1/4"	80-982-17-04
0.3 mm	20 cm/8"	80-982-20-04
<b>blunt</b>		
0,6 mm	12 cm/4 3/4"	80-984-12-04
0,6 mm	17 cm/6 1/4"	80-984-17-04
0,6 mm	20 cm/8"	80-984-20-04
<b>blunt</b>		
1.0 mm	12 cm/4 3/4"	80-986-12-04
1.0 mm	17 cm/6 1/4"	80-986-17-04
1.0 mm	20 cm/8"	80-986-20-04
1.0 mm	23 cm/9"	80-986-23-04
<b>blunt</b>		
2.0 mm	17 cm/6 1/4"	80-988-17-04
2.0 mm	20 cm/8"	80-988-20-04
2.0 mm	23 cm/9"	80-988-23-04
2.0 mm	30 cm/11 3/4"	80-988-30-04



revolutionary  
ergonomic  
design

Bipolar forceps, KLS Martin nonStick red, angled



Tip width	Length	Item No.
<b>pointed</b>		
0.3 mm	12 cm/4 3/4"	80-983-12-04
0.3 mm	17 cm/6 1/4"	80-983-17-04
<b>blunt</b>		
0,6 mm	12 cm/4 3/4"	80-985-12-04
0,6 mm	17 cm/6 1/4"	80-985-17-04
0,6 mm	20 cm/8"	80-985-20-04
<b>blunt</b>		
1.0 mm	12 cm/4 3/4"	80-987-12-04
1.0 mm	17 cm/6 1/4"	80-987-17-04
1.0 mm	20 cm/8"	80-987-20-04
1.0 mm	23 cm/9"	80-987-23-04
<b>blunt</b>		
2.0 mm	17 cm/6 1/4"	80-989-17-04
2.0 mm	20 cm/8"	80-989-20-04
2.0 mm	23 cm/9"	80-989-23-04
2.0 mm	25 cm/9 3/4"	80-989-25-04
2.0 mm	28 cm/11"	80-989-28-04
2.0 mm	30 cm/11 3/4"	80-989-30-04

1/1

80-983-12-04  
12 cm/4 3/4"  
pointed, 0.3 mm

Example illustration

revolutionary  
ergonomic  
design




**Bipolar forceps, KLS Martin nonStick red, bayonet-shaped**

Tip width	Length	Item No.
<b>pointed</b>		
0.3 mm	17 cm/6 1/4"	80-990-17-04
0.3 mm	20 cm/8"	80-990-20-04
0.3 mm	23 cm/9"	80-990-23-04
0.3 mm	25 cm/9 3/4"	80-990-25-04
<b>blunt</b>		
0,6 mm	17 cm/6 1/4"	80-991-17-04
0,6 mm	20 cm/8"	80-991-20-04
0,6 mm	23 cm/9"	80-991-23-04
0,6 mm	25 cm/9 3/4"	80-991-25-04
<b>blunt</b>		
1.0 mm	17 cm/6 1/4"	80-992-17-04
1.0 mm	20 cm/8"	80-992-20-04
1.0 mm	23 cm/9"	80-992-23-04
1.0 mm	25 cm/9 3/4"	80-992-25-04
<b>blunt</b>		
2.0 mm	17 cm/6 1/4"	80-993-17-04
2.0 mm	20 cm/8"	80-993-20-04
2.0 mm	23 cm/9"	80-993-23-04
2.0 mm	25 cm/9 3/4"	80-993-25-04

Example illustration



Example illustration




upward-bent

downward-bent

revolutionary  
ergonomic  
design

**Bipolar forceps, KLS Martin nonStick red, bayonet-shaped, angled**



Tip width	Length	Item No.
<b>blunt</b>		
0.6 mm	23 cm/9"	80-997-23-04
0.6 mm	25 cm/9 3/4"	80-997-25-04
<b>blunt</b>		
0.6 mm	23 cm/9"	80-996-23-04
0.6 mm	25 cm/9 3/4"	80-996-25-04
<b>blunt</b>		
1.0 mm	20 cm/8"	80-994-20-04
1.0 mm	23 cm/9"	80-994-23-04
1.0 mm	25 cm/9 3/4"	80-994-25-04
<b>blunt</b>		
1.0 mm	20 cm/8"	80-995-20-04
1.0 mm	23 cm/9"	80-995-23-04

Bipolar forceps,  
KLS Martin  
nonStick red,  
angled



$\frac{1}{1}$   
80-998-23-04  
23 cm / 9"

$\frac{1}{1}$   
80-999-23-04  
23 cm / 9"


$\frac{1}{1}$   
80-999-24-04  
23 cm / 9"

$\frac{1}{1}$   
Meuser  
80-998-24-04  
18 cm / 7"

**Connecting cables for bipolar forceps  
with angled plugs**




80-291-40-04  
4 m/13 ft.  
Connecting cable for bipolar instruments  
for KLS Martin and Berchtold HF units



80-286-40-04  
4 m/13 ft.  
Connecting cable for bipolar instruments  
for maxium® "e"-version / Erbe HF units ICC  
and ACC/VIO



80-293-40-04  
4 m/13 ft.  
Connecting cable for bipolar instruments  
for maxium® "i"-version, ME MB2 "i"-version /  
Valleylab HF units




## Storage and transport

To store and transport your nonStick red forceps safely and securely, we recommend using the new KLS Martin mesh tray system.

A lid is optionally available for each tray size, to be easily locked in place by simply folding down the tray handles.

The forceps can be fixed in place in the tray fully in accordance with your needs.



Fixing elements are available as 6-item strips with three holding clips. However, you can easily cut up the strips with a pair of scissors – just as you need it.



For forceps with a maximum length of 23 cm, we recommend the following configuration:

Item No.	Designation	Quantity
55-804-25-01	Mesh tray, 243 x 255 x 53 mm	1
55-805-28-01	Lid, 243 x 255 mm	1
<b>Fixing elements</b>		
55-806-02-04	Clip for fixing element (10/pack)	1
55-806-07-04	Fixing element (6/pack)	2

For forceps longer than 23 cm, a ¾-size tray should be used.

Item No.	Designation	Quantity
55-804-42-01	Mesh tray, 410 x 255 x 53 mm	1
55-805-45-01	Lid, 410 x 255 mm	1
<b>Fixing elements</b>		
55-806-02-04	Clip for fixing element (10/pack)	1
55-806-07-04	Fixing element (6/pack)	2



Item No.	Designation	Quantity
55-804-19-01	marTray®, mini for Arthroscopy 277 x 84 x 36 mm	1



For an overview of the whole range of KLS Martin HF products, please refer to our Catalog "Accessories for electrosurgery".



The respectful treatment of people, society and the environment has always been a matter of course to us. By design, most of our products can be recycled and used over many years, thus saving resources and reducing waste. Environment-friendly and recyclable materials are used also for our production processes, and energy and water consumption monitoring is a top priority as well. A case in point: our largest manufacturing site includes a heat recovery plant for efficient energy utilization. These and many other measures are clear evidence of our environmental commitment.

For further information, please refer to separate GoGreen brochure.

## KLS Martin Group

### KLS Martin Australia Pty Ltd.

Sydney · Australia  
Tel. +61 2 9439 5316  
australia@klsmartin.com

### KLS Martin do Brasil Ltda.

São Paulo · Brazil  
Tel. +55 11 3554 2299  
brazil@klsmartin.com

### KLS Martin Medical (Shanghai) International Trading Co., Ltd

Shanghai · China  
Tel. +86 21 5820 6251  
info@klsmartin.com

### KLS Martin India Pvt Ltd.

Chennai · India  
Tel. +91 44 66 442 300  
india@klsmartin.com

### KLS Martin Italia S.r.l.

Milan · Italy  
Tel. +39 039 605 67 31  
info@klsmartin.com

### KLS Martin Japan K.K.

Tokyo · Japan  
Tel. +81 3 3814 1431  
info@klsmartin.com

### KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia  
Tel. +604 261 7060  
malaysia@klsmartin.com

### KLS Martin de México, S.A. de C.V.

Mexico City · Mexico  
Tel. +52 55 7572 0944  
mexico@klsmartin.com

### KLS Martin Nederland B.V.

Huizen · Netherlands  
Tel. +31 35 523 45 38  
infonl@klsmartin.com

### KLS Martin SE & Co. KG

Moscow · Russia  
Tel. +7 499 792 76 19  
russia@klsmartin.com

### KLS Martin Taiwan Ltd.

Taipei · Taiwan  
Tel. +886 2 2325 3169  
taiwan@klsmartin.com

### KLS Martin SE & Co. KG

Dubai · United Arab Emirates  
Tel. +971 4 454 16 55  
middleeast@klsmartin.com

### KLS Martin UK Ltd.

Reading · United Kingdom  
Tel. +44 118 467 1500  
info.uk@klsmartin.com

### KLS Martin LP

Jacksonville · Florida, USA  
Tel. +1 904 641 77 46  
usa@klsmartin.com

### KLS Martin SE Asia Sdn. Bhd.

Hanoi · Vietnam  
Tel. +49 7461 706-0  
info@klsmartin.com

Do you know how to get all important information  
about the KLS Martin Energy Devices?

Please download KLS Martin App Energy Devices!  
The App is available for Android and iOS.



<https://itunes.apple.com/de/app/klsmartin-energy-devices/id1198171415?l=de&ls=1&mt=8>



<https://play.google.com/store/apps/details?id=com.klsmartin.energydevices>

### KLS Martin SE & Co. KG

#### A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany  
PO Box 60 · 78501 Tuttlingen · Germany  
Tel. +49 7461 706-0 · Fax +49 7461 706-193  
info@klsmartin.com · www.klsmartin.com