



FlexPen 2.0

Ergonomic and flexible!


Fiber Holder **FlexPen** 2.0 Ergonomic and flexible!

The new fiber holder FlexPen 2.0 from KLS Martin combines ergonomics, easy handling and flexible use in one instrument.

FlexPen 2.0 is a modular system for the guidance and fixation of a laser fiber or a waveguide and is comprised of a handle, a Tuohy-Borst adapter and various attachments with bending mandrel.


Maximum ergonomics

Due to its grip-friendly design and ribs the instrument is extremely ergonomic to hold optimally so it allows comfortable work.



Many different functions

The Luer lock attachment makes it possible to connect irrigation or suction tubes.








Easy handling and reprocessing

The Tuohy-Borst adapter, which fixates the fiber securely and firmly, is disposed of after use – so there are no more small, barely accessible parts that can get lost during disassembly or reprocessing!



Flexible use options:

Not only do the different attachments make it possible to use the instrument in various disciplines and indications – all the available fiber diameters are easy to insert and fixate.

	79-311-00-04 Fiber holder handle FlexPen 2.0
	79-311-03-04 Tuohy-Borst adapter (pack of 25, sterile, single use)
	79-311-02-04 7-cm attachment tip with suction (2-piece: suction tube and internal fiber guide tube)
	79-311-10-04 10-cm attachment tip (incl. bending mandrel)
	79-311-30-04 30-cm attachment tip for waveguide (incl. bending mandrel)

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co., Ltd

Shanghai · China
Tel. +86 21 5820 6251
info@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

KLS Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
info@klsmartin.com

KLS Martin Japan K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
info@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 261 7060
malaysia@klsmartin.com

KLS Martin de México, S.A. de C.V.

Mexico City · Mexico
Tel. +52 55 7572 0944
mexico@klsmartin.com

KLS Martin Nederland B.V.

Huizen · Netherlands
Tel. +31 35 523 45 38
infonl@klsmartin.com

KLS Martin SE & Co. KG

Moscow · Russia
Tel. +7 499 792 76 19
russia@klsmartin.com

KLS Martin Taiwan Ltd.

Taipei · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

KLS Martin SE & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

Reading · United Kingdom
Tel. +44 118 467 1500
info.uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Hanoi · Vietnam
Tel. +49 7461 706-0
info@klsmartin.com

Do you know how to get all important information
about the KLS Martin Energy Devices?

Please download KLS Martin App Energy Devices!
The App is available for Android and iOS.



<https://itunes.apple.com/de/app/klsmartin-energy-devices/id1198171415?l=de&ls=1&mt=8>



<https://play.google.com/store/apps/details?id=com.klsmartin.energydevices>

KLS Martin SE & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
PO Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com